



Recombinant Mouse Zinc finger C4H2 domain-containing protein (Zc4h2)

Product Code	CSB-EP734664MO-B
Abbreviation	Zc4h2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68FG0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MADEQEIMCK LESIKEIRNK TLQMEKIKAR LKAEFEALES EERHLKEYKQ EMDLLLQEKM AHVEELRLIH ADINVMENTI KQSENDLNKL LESTRRLHDE YKPLKEHVDA LRMTLGLQRL PDLCEEEEEKL SLDYFEKQKA EWQTEPQEPP IPESLAAAAA AAQQLQVARK QDTRQTATFR QQPPPMKACL SCHQQIHRNA PICPLCKAKS RSRNPKKPKR KQDE
Source	E.coli
Target Names	Zc4h2
Protein Names	Recommended name: Zinc finger C4H2 domain-containing protein Alternative name(s): Hepatocellular carcinoma-associated antigen 127 homolog
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.