



# Recombinant Rat Conserved oligomeric Golgi complex subunit 6 (Cog6)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP734667RA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Abbreviation</b>      | Cog6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Uniprot No.</b>       | Q68FP9                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Sequence</b>          | MADTSGEVAA VPASGAANGL SNGAGATPAQ PNNPLSRKLH KILETRLND<br>KEMLEALKAL SAFFVENSLR TRRNLRGDIE RRSLAINEEF VSIFKDVKEE<br>LESINEDVQA MSSCCQDMTS RLQAAKEQTQ DLIVKTTKLQ AENQRLEIRA<br>QVADAFLSKF QLTSDMTLL RGTRGGPVTE DFFKALGRVK QIHNDVKVLL<br>RTNQQTAGLE IMEQMALLQE TAYERLYRWA QSECRALTQE SCDVSAVLTQ<br>AMEALQDRPV LYKYTLDEFG TARRSTVVRG FIDALTRGGP GGTPRPIEMH<br>SHDPLRYVGD MLAWLHQATA SEKEHLEALL KHVTAQGVKE NIQEVVGHIT<br>EGVCRPLKVR IEQVILAEPG AVLLYKISNL LKFYHHTISG IVGNSAATLL<br>TTIEEMHLLS KKIFFSSLSL HANKLMDKIE LPPDGLGPSS ALSQTLTLLR<br>DVLASHDSSV VPLDARQADF VQVLSCVLDL LLQMCTVSAS NLGTADMATF<br>MVNSLYMMKT TLALFEFTDR RLEMLQFQIE AHLDTLINEQ ASYVLTRVGL<br>SYIYNTIQQH RPDQGSLSM PNLDVSLKA AMAQFDRYLS APDHLLMPQL<br>NSLLSATVKE QIIKQSTELV CRAYGEVHAA VMNPVNAYKD PESILHRSPE<br>QVKTLLS |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Target Names</b>      | Cog6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Names</b>     | Recommended name: Conserved oligomeric Golgi complex subunit 6 Short name= COG complex subunit 6 Alternative name(s): Component of oligomeric Golgi complex 6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Expression Region</b> | 1-657                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.