



Recombinant Pygathrix roxellana Insulin-like growth factor I (IGF1)

Product Code	CSB-YP734672PGAL
Abbreviation	IGF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68LC0
Product Type	Recombinant Protein
Immunogen Species	Rhinopithecus roxellana (Golden snub-nosed monkey) (Pygathrix roxellana)
Purity	≥85% (SDS-PAGE)
Sequence	GP ETLCGAELVD ALQFVCGDRG FYFNKPTGYG SSSRRAPQTG IVDECCFRSC DLRRLEMYCA PLKPAKSA
Source	Yeast
Target Names	IGF1
Protein Names	Recommended name: Insulin-like growth factor I Short name= IGF-I Alternative name(s): Somatomedin
Expression Region	49-118
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.