



Recombinant *Oryza sativa* subsp. *japonica* Tubby-like F-box protein 9 (TULP9)

Product Code	CSB-MP734731OFG
Abbreviation	TULP9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68Y48
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MALWRCSSSW LSSVSRSSGG VGGGESKVSP EIAPVSGGEG EGEEEEEEEE RWSRLLPELL TEIMRRVDAG AERWPPRRDV VACACVCRRW RDAAVSVVRP PLECGRITFP SSLKQPGPRD APMHCFIRRN KKNSTFYLYL SLTQALTDKG KFLLAARRFR NGAHTEYIIS YDCDDLFPGS NSYVGKLRSD FLGTFKFIYD SQPPYDGAKP SRSQSSRRFA SKQINPNVSG GNYEVGQVSY KFNFLKSRGP RRMQCNIQCP VGQSTASDPL KKLISTSSPL ALRNKAPRWH EHLQCWCLNF HGRVTVASVK NFQLVAPAGT SDPWGIADEE TVILQFGKIE DDAFTMDYRQ PLSAFQAFAI CLTSFGTKLA CE
Source	Mammalian cell
Target Names	TULP9
Protein Names	Recommended name: Tubby-like F-box protein 9 Short name= OsTLP9 Alternative name(s): Tubby-like F-box protein 10 Short name= OsTLP10
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.