



Recombinant Propionibacterium acnes DNA polymerase IV (dinB)

Product Code	CSB-BP734782PSG
Abbreviation	dinB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6A763
Product Type	Recombinant Protein
Immunogen Species	Propionibacterium acnes (strain KPA171202 / DSM 16379)
Purity	≥85% (SDS-PAGE)
Sequence	MRSTASILHL DMDAFFASVE QRDKPSLVGK AVIVGGVGGR GVVATASYEA RKFGVHSAMA GSQARRLAPN AAFLSGRFES YRESSKVVMA TLRELSPLVE PLSLDEAFVD LEAGDVDVEN LDALQQVAVE LRARVAERTE GLSCSVGIGS SKFMAKVASE MAKPQRGHTD QIALIAPGTE ADTIAPLSAR AIPGVGPVTA ERLEKLGLRT IADVRAARES ELVHELGHAS GASLVALSRA RDDRPPVASR IAKSISVEDT FEHDLTHREE CRQIVERHAG LVCGRLLKSG QFARTVTLKA KMAFTVWSR SVTLTGATDS PERISQLAVR MLDSLRLREG VRLGVLGVSN FATSAQEELF VVDDEGRLLD EPEVTEEITP IQPGPFGRRR SFDGRNWPPG VDLIHDEYGR GWVWGSGHGI VTCRFEYRGS EIGRVKSVQP DDEALHPAHP LPLAWQPPQE DVPSE
Source	Baculovirus
Target Names	dinB
Protein Names	Recommended name: DNA polymerase IV Short name= Pol IV EC= 2.7.7.7
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.