



Recombinant Leifsonia xyli subsp. xyli Diaminopimelate epimerase (dapF)

Product Code	CSB-YP734814LAAC
Abbreviation	dapF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AE05
Product Type	Recombinant Protein
Immunogen Species	Leifsonia xyli subsp. xyli (strain CTCB07)
Purity	≥85% (SDS-PAGE)
Sequence	MATLHFTKGQ GTGNDFVLYA DPDGTLPLSA GQIAEICDRH FGVGADGVIR AVRSKRLPEG AEALADDEAA EWFMDYYNAD GSAAEMCGNG IRVYVRYLIE SRLVQLADGD TLPIGTRSGV RDVQRTLSGF QVDLGRWRLA GGEPLVRAKE LPVARLGLGI DLGNPHVVVA LADESELESA DLTYIPHLEP APQGGANVEF VVPHEPLVKD GVGRIRMRVH ERGSGETLSC GTGAAAAALA VRHWAGENAP GQWWVDVPGG TVGVRMFPTD DGEHVALSGP AELVYTGALE LA
Source	Yeast
Target Names	dapF
Protein Names	Recommended name: Diaminopimelate epimerase Short name= DAP epimerase EC= 5.1.1.7
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.