



Recombinant Rat Pre-mRNA-splicing factor 38B (Prpf38b)

Product Code	CSB-YP734892RA
Abbreviation	Prpf38b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AXY7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	<p>ANNSPALTG NSQPQHAAA AVVQQQQQCG GGGATKPAVS GKQGNVLPW GNEKTMNLNP MILTNILSSP YFKVQLYELK TYHEVVDEIY FKVTHVEPWE KGSRKTAGQT GMCGGVRGVG TGGIVSTAFCLLYKLFTLKL TRKQVMGLIT HTDSPYIRAL GFMYIRYTQP PTDLWDWFES FLDDEEDLDV KAGGGCVMTI GEMLRSLTK LEWFSTLFPRI PVPVQKNID QQIKTRPRKI KKDGGKEGVVE IDRHIERRRS RSPRRSLSPR RSPRRSRRSRSHHRDGHGSSSS FDRELEREKE RQRLEREAKE REKERRRSRS LDRGLDRRRSRSRERHRSRS RSRDRKGDRR DRDREREKEN ERGRRRDRDY DKDRGNDRER DRERSRERSK EQRSRGESEE KKHKEDKEDR RHRDDKKEK KKHSRSRSRE RKHSRSRSR NAGKRSRSRS KDKASKHKNE SKEKSNKRSR SGSQGRGTSV EKRKREHSPS REKSRKRSRS QDRSHKRDHD SKDQSDRQDH QRSQSVEPES QEKEHKNKDE TV</p>
Source	Yeast
Target Names	Prpf38b
Protein Names	Recommended name: Pre-mRNA-splicing factor 38B
Expression Region	2-542
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.