



Recombinant Rat Protein FAM203A (Fam203a)

Product Code	CSB-MP734899RA
Abbreviation	Fam203a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AY79
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	RGPGVGSGL SGERGLSGKE AGADTEVAEL LPFLVPGARA DLQAAAAQHV LALTGVGSGR TLLAGQTALL RALIDLAVAP IPAPSRDASR ALVNLAADPD LHWQMLAADP ELPVRLRCV LDPQWPWAE LAAVLANLSR EPAPCAVLTE KLVAAPERL GLERLVSALC MPSYNASAPL HYLGPLLSNL SQQAEVRAFL LDPDRCVVQR LLPLTQYTDS SVRRGGVVG TLRNCCFEHRH HKWLLGPQVD ILPFLLLPLA GPEEFSEEEM EELPVDLQYL SPDKQREPDA DIRKMLIEAI MLLTATAPGR KQVRDQGAYL ILRELYAWES EPDVRMACEK LIQVLIGDEP EVGMENLLEV QVPEDVERQL QELDQQEQQE LAQERRDKGA PRT
Source	Mammalian cell
Target Names	Hgh1
Protein Names	Recommended name: Protein FAM203A Alternative name(s): Brain protein 16
Expression Region	2-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.