



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized protein YMR135W-A (YMR135W-A)

Product Code	CSB-MP734918SVG
Abbreviation	YMR135W-A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6B0X1
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MAEYLATCSF HSSFVKHFFP YVKSFPPLP WPLQDPLFVV SLLRVDIVIS ICYHTQIYLH PADICLYCPF ACNLMTKLHM LASRKMMYHQ NVSCNEPGHR AGRVRKARST LIVINSNTME RLPFTRNGSG QQSNKLRDPK KGRTHKPKPS EKHKKNKTGK KGAQEKTHRS RSSRKGN
Source	Mammalian cell
Target Names	YMR135W-A
Protein Names	Recommended name: Putative uncharacterized protein YMR135W-A
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.