



Recombinant *Debaryomyces hansenii* RuvB-like helicase 1 (RBV1)

Product Code	CSB-YP734950DIS
Abbreviation	RBV1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BI60
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MVQINEVKDT QTRESRTAAH THIKGLGLDE HGIAKRVEGG FVGQSDAREA CGIIVDLIKS KRMSGKAILL AGAPGTGKTA LALAIQELG PKVPCPIVG SELFSAEIKK TAALMENFRR AIGLRIKETK EVYEGEVIEL TPEEAENPLG GYGKTISHVI VGLKTAKGTK SLRLDPSIYE SIQKERVSVG DVIIYEANTG SVKRVGRSDA YATEFDLEAE EYVPLPKGEV HKKKEIVQDV TLHDLDVANA RPQGGQDVLS MMGQLLKPRK TEITDKLRSE VNKVVSKYID QGVAELIPGV LFIDEVNMLD IECFTYLNRA LESSIPIVV LASNRGMTTI RGTDDDKKSP HGCPADLIDR LLIVRTL PYN QEEIKIIISK RATLENLIVT PDALDKLSLH GINNSLRYAL QLLAPAGVLS KTAGRNEITS EDIEECEILF LDSRRSIKIL EETKGYL
Source	Yeast
Target Names	RBV1
Protein Names	Recommended name: RuvB-like helicase 1 EC= 3.6.4.12
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.