



Recombinant Debaryomyces hansenii Pre-rRNA-processing protein IPI1 (IPI1)

Product Code	CSB-BP734965DIS
Abbreviation	IPI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BLC4
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MGSKRKKA EK QKDFVKAKLR VGKTA AKPDN HTDTSFIAKS ISIPNQ TINK KTSNTEKKYE VDLAHHLSLT KHHSDVTRKE VLNYIEQHLP SNPSLYKDIL TSIVPLIIDQ SQNVRNALTS LLSACATQQV GLLDLHIRSI ILFIHSAMSH IKPDIRNSST KFLSVLIDHA TESLVRSYFI KTLKSYFTLL SWTLTNDKKA VSLAITTSSS IGGPSKKARI HHLSILRAFL SAALFPISSN ERKLDYSKIK MIHPESYKYL LASSTQPFAS LKLFVQEV PK QKNTDIQQNT KKDDNAFSLN DLDTVSTEDI DTRRKVMLDV FMAPMLRNLK NLIKEGGEVG REAHSCMKVL EQLQSETK
Source	Baculovirus
Target Names	IPI1
Protein Names	Recommended name: Pre-rRNA-processing protein IPI1
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.