



Recombinant Debaryomyces hansenii ATP-dependent rRNA helicase RRP3 (RRP3)

Product Code	CSB-MP734982DIS
Abbreviation	RRP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BQ61
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MSVKVDGMIN KSKTHSKKL DARALAEKIK KNALKQQAQE EPAKEENVAE NFDTVVDPTA ELKFKSFNEL KLIPELLEAI QQMKFTKPTP IQSEAIPHAL EGKDIIGLAQ TGSGKTAFA IPILQALWEA QAAYYGLVLA PTRELAYQIK ETFDALGSSM GLRSVCIVGG MDMMDQARDL MRKPHILVAT PGRIMDHLEH TKGFSLKNLK YLVMDEADRL LDMDFGPALD KILKIPTQR TTYLFSATMT NKIAKLQRAS LHPVVRVAVS NKYQTADNLV QSMMLVSDGY KNTYLIHLLN EFLGKSIIF TRTCAHSQRT ALLARILGFS AVPLHGQLTQ AQRLGSLNKF KAGKANILIA TDVAARGLDI PSVDIVINYD IPTDSKAYIH RVGRTARAGK SGKSISLITQ YDLEMYLRIE SVLGKKLPKD PSPSKAVLDT LHVHVDRASA EAIRQTKDFH EKRNPCKNRD DRDREER
Source	Mammalian cell
Target Names	RRP3
Protein Names	Recommended name: ATP-dependent rRNA helicase RRP3 EC= 3.6.4.-
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.