



Recombinant *Yarrowia lipolytica* Autophagy-related protein 3 (ATG3)

Product Code	CSB-BP735039YAP
Abbreviation	ATG3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C4Q9
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	>85% (SDS-PAGE)
Sequence	MEAHKGGKQDI FIVQVLITPT FVELNHIYPP IRLLSPPSIS SNPSIKFDKT QAPPTLFIYI ASYNTPDYQT TLHKSTMLHV RAAASSLREY LTPVSNTSTF RTTGEITPDE FVKAGDYLVE KFPTWSWASA SKSKVRDFLP PDKQQLVTRH VPSHVRASTV SGPVTLGEEE DGWTSFGVAN TKDADGDDTE EIAEIASDSD EELDDDDDDA AEAPATATDH RTYNLYIAYS TSYRVPKMFL SGYSPEGSPL TPEDMFEDII PEYRDKVTI ERPTFQDNIT MVAIHCKHA NVMRVLMERV EAKGDKDITR GVAKLGVADA DDGEGEEWE EVENSAMRVD QYLVTFKFI ASVTPGIEHD YTMSAL
Source	Baculovirus
Target Names	ATG3
Protein Names	Recommended name: Autophagy-related protein 3 Alternative name(s): Autophagy-related E2-like conjugation enzyme ATG3
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.