



# Recombinant *Yarrowia lipolytica* Chromatin structure-remodeling complex subunit SFH1 (SFH1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP735058YAP-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Abbreviation</b>      | SFH1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                               |
| <b>Uniprot No.</b>       | Q6C9N2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | <i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) ( <i>Candida lipolytica</i> )                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Sequence</b>          | MGEKPNRIGR ANISYINHDH IFLPSFQPSY PSRQLHPLIS SQNKVHSSIV<br>HLPKVASPSS PETQQHSQAP LFITHIMSGP PQALASGFAA RVREGGTTLY<br>INTTPMLRSA RHNTAVNYAE FEEDFDANDF EDDDDDDQSQ RESRDGSEEA<br>EGDEDGKKE EQDKFAGLKA PLVSNEPKRA APPVRPVMYP QEVLEELSQV<br>KEPTLIPIRV AVENIDVFRV QDFFLWDADE KILTPEQFAT LTCADLDVPI<br>GYSAQMSAQI KKQLAEYTAA PALPKDVEVH VIVELAVTV D KIVYEDKFEW<br>DLSGEYATPQ EFARTVVQDL GLGQEFYPAI TYQLYETLGK LQKAWLERSI<br>PLDVDNRAAF GLEAGLRVDQ DNLGESWVPR VEEMTPEEMQ<br>KREMERDRSS RRLKRESARM AEVPYVDLDS LYSRKRRRRF<br>DEDSRSGSPM W |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>      | SFH1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: Chromatin structure-remodeling complex subunit SFH1<br>Alternative name(s): RSC complex subunit SFH1 SNF5 homolog 1                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Expression Region</b> | 1-441                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                      |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.