



Recombinant *Yarrowia lipolytica* Chromatin structure-remodeling complex subunit SFH1 (SFH1)

Product Code	CSB-MP735058YAP
Abbreviation	SFH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C9N2
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MGEKPNRIGR ANISYINHDH IFLPSFQPSY PSRQLHPLIS SQNKVHSSIV HLPKVASPSS PETQQHSQAP LFITHIMSGP PQALASGFAA RVREGGTTLY INTTPMLRSA RHNTAVNYAE FEEDFDANDF EDDDDDDQSQ RESRDGSEEA EGDEDGKKE EQDKFAGLKA PLVSNEPKRA APPVRPVMYP QEVLEELSQV KEPTLIPIRV AVENIDVFRV QDFFLWDADE KILTPEQFAT LTCADLDVPI GYSAQMSAQI KKQLAEYTAA PALPKDVEVH VIVELAVTVD KIVYEDKFEW DLSGEYATPQ EFARTVVQDL GLGQEFYPAI TYQLYETLGK LQKAWLERSI PLDVDNRAAF GLEAGLRVDQ DNLGESWVPR VEEMTPEEMQ KREMERDRSS RRLKRESARM AEVPYVDLDS LYSRKRRRRF DEDSRSGSPM W
Source	Mammalian cell
Target Names	SFH1
Protein Names	Recommended name: Chromatin structure-remodeling complex subunit SFH1 Alternative name(s): RSC complex subunit SFH1 SNF5 homolog 1
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.