



Recombinant Yarrowia lipolytica Pre-mRNA-splicing factor SLU7 (SLU7)

Product Code	CSB-EP735059YAP-B
Abbreviation	SLU7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C9T0
Product Type	Recombinant Protein
Immunogen Species	Yarrowia lipolytica (strain CLIB 122 / E 150) (Yeast) (Candida lipolytica)
Purity	≥85% (SDS-PAGE)
Sequence	MSKPATTLDG LQDGSVELQK RLSAKRDAEG ELKHEEDRKL RKEGKTQTPV TNQEGDVLNP YIPRFMSSAP WYAKAGDDDP SASSLDHLKK DKTDSKKEEW YARGQDRSSL ERPKKFRKGA CENCGAMTHK KKDCMERPRK VGA KYKADNL QADDVISNPH MSWDSKRDRW NGYDAEDFKR VIDEHNMTTE LQKKLHGDAP AETDAEKAAA VEEAIRKSTR TLRLREDTAV YLKDLNSETV YNPGSRTFRT ADEGYIDKDG MFVRHTTGEG KEMEELTRIA EAESALGNQI HLHANPTAAA IAAKRIKEEA EKAKDAEKRA LEEKYGKQDH VRVRPKEYTS VVTAKKPTVI RAKNAQGI AV SKYEEDVFPG NHTSVWGSYY DKETAKWGYA CCHCLIKNGY CVGEKGK KKK GAMDTV
Source	E.coli
Target Names	SLU7
Protein Names	Recommended name: Pre-mRNA-splicing factor SLU7
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.