



Recombinant *Yarrowia lipolytica* Survival factor 1 (SVF1)

Product Code	CSB-YP735066YAP
Abbreviation	SVF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CBN5
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MLKWWQGGVS ALTGSAEPEY GAEAFESVSK SVDGKNPYGS LTLPELNWLK PAGSNVETQT FYFFSKDGKQ FGHVQIIHSN PLGIQNTAQF TFRLSNADKP EENVWSSTNV EEFEVGPVF TAKGPMGVDN IKLVPSVDDD GKTTEAFHVK IHVSDEVKLD LTFTRSVGEF KIGPTGQSKY GEDLNNPWGA MRHSFWPRAS VSGTIEVANK DKMELDNDFG FFVHATQGMK PHHAASKWNF ITYQSPKYSA VVM DYTTTPS YGSSRV TIGG VASAEKLLFT SANATVEHVE TKEDEEVGWQ VPTLIKFECE GPKVDVADED VDKPENRGKA VLQGNLTHLE RMDVMAEIPA FLKSLASTVS GTKPYIYQFR NEMELELEID GEKIQDKGQI FCEATFIS
Source	Yeast
Target Names	SVF1
Protein Names	Recommended name: Survival factor 1
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.