



Recombinant Kluyveromyces lactis Protein DML1 (DML1)

Product Code	CSB-BP735148KBK
Abbreviation	DML1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CU61
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MREIINVSVS HRSNHLITQF YNCLEPLLHD ADQENDVFLN PNIDKVSKTV SYTPRALLWD AKLGNGSLGT YQYVSENDYA DTLTSEQGGAT AKTAHRVQTH DRIRKSPYQL ALDQGATVLP KINDEIAKYW SDYSKLIYDP SSFNTLQDWY HDAANQQKAP NFQNLQRQVYF DNYETGSNQF RENYSNEFFD SNLHQQLEKC DSLQGFNIIT ELDNGWGGFS SSMLLELKDE LPKVSYHTYG WNQDDVCSLK EPVHSTKTKF QMLCNKIRAT IALSQESDLF FPLYANTKEA FWRSTGETVL LFDSVNSVFH GSRKQAATNS NMRKMSHLTN ALTLGESKRN IVSKLNVDEY NSFSFYDRMP LYKNSKKPSH TFTQCEIDRV NHKESSMFHF TTYPWTNSDT IPDTHHSCEAE SVIGVTEKPR DVLKWTWEDLV SRYFRYDSDR EEIKDHLGTL SLEYEHGWYD DDDSDLDL
Source	Baculovirus
Target Names	DML1
Protein Names	Recommended name: Protein DML1
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.