



Recombinant Kluyveromyces lactis Nascent polypeptide-associated complex subunit alpha (EGD2)

Product Code	CSB-BP735161KBK
Abbreviation	EGD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CWG6
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MSEIAQDSNV SIFSKNEKKA RELILKLGLK PIQGVSRVTF KKKDGGQIFAI DRPEVYKSQG GNFFVFGGEAK VDDFTQRLAK AQESLQQNEG VLPAGQDAVS KDPQSIQADM QAAADSATDK PSADDAVDDA EVDETGLNAD DIELVMQQAG VPRAKAAKAL KEHDSDIVNA IMALSG
Source	Baculovirus
Target Names	EGD2
Protein Names	Recommended name: Nascent polypeptide-associated complex subunit alpha Short name= NAC-alpha Alternative name(s): Alpha-NAC
Expression Region	1-176
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.