



Recombinant *Erwinia carotovora* subsp. atroseptica 3-isopropylmalate dehydratase small subunit (leuD)

Product Code	CSB-BP735194EAAB
Abbreviation	leuD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D0G5
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MAKFTQHTGL VVPLDAANVD TDAIPKQFL QKVTRTGFGQ HLFHDWRFLD DAGQQPNPEF VLNKPHYKGA SILLARENFG CGSSREHAPW ALTDYGFSSV IAPSFADIFY GNSFNNQLLP VKLSDEEVDE LFKLVDEQEG ITFTVDLENQ VVQADSKRYP FEIDSFRRHC MINGLDSIGL TLQHEASITE YEKNQPAFLN
Source	Baculovirus
Target Names	leuD
Protein Names	Recommended name: 3-isopropylmalate dehydratase small subunit EC=4.2.1.33 Alternative name(s): Alpha-IPM isomerase Short name= IPMI Isopropylmalate isomerase
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.