



# Recombinant *Erwinia carotovora* subsp. *atroseptica* Elongation factor Ts (tsf)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP735246EAAB                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Abbreviation</b>      | tsf                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q6D8E2                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) ( <i>Erwinia carotovora</i> subsp. <i>atroseptica</i> )                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MAEITASLVK ELRERTAAGM MECKKALVEA NGDIELAIEN MRKSGAIKAA<br>KKAGNVAADG VIKTKIDGNY AVILEVNCQT DFVAKDGGFQ AFADKVLDA<br>VAGKITDVDV LKAQFEEERV ALVAKIGENI NIRRVSALLEG EVLGNYQHGA<br>RIGVLVAAKG ADEELVKHLA MHVAASKPEF VKPEDVSAEV VAKEYQVQLE<br>IAMQSGKPKE IAEKMVEGRM KKFTGEVSLT GQPFVMDPAK SVGQLLKEHN<br>ADVTNFIRFE VGEGIEKVET DFAAEVAAMS KQS                                                                      |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Target Names</b>      | tsf                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Names</b>     | Recommended name: Elongation factor Ts Short name= EF-Ts                                                                                                                                                                                                                                                                                                                                                    |
| <b>Expression Region</b> | 1-283                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |