



Recombinant Danio rerio Ubiquitin-conjugating enzyme E2 variant 3 (uevld)

Product Code	CSB-EP735269DIL
Abbreviation	uevld
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DBY5
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDVNSEPVKK VLSKYKFRDV AVEELQKVYR VYPDMKIMAG TYTSSDSLQK DLLKLVGNIP VVYQGRSYNI PILLWLLDSF PFTPPICYLR PTSSMVIREG KHVDSKGRIH LPALHNWDHP KSSVNALLAE MIGKFEEPP LGTKSSAHGD TPSNLLDYVS NLTITEGGNR PDQEVKSVI GGGDLGIAAV LSIMAKSCVD KLVLIDIPEN STKGGTMDLE IFSLPKVEVS KDLSASAGSK VLVITANAWS DEQSYLSVVQ TNVDMYRGII PRLAQLSPNA VLVIASQPVD VMTHVAWRQS HLLPTQVIGV GCNLD SQRLS HIINISLVAN NTGKQAWVIG ELSNKVAVW GNMGP GTDQL QALTPVSNST KPLMDRAFEM IKGRGQRSWS VGLSIADITH SIVTNQKKVH SVTTLAEGWG GIGSKVFLSL PCILGETGST RLPGVALGSE DEMKLRESVA CQVNLFMQLR L
Source	E.coli
Target Names	uevld
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 variant 3 Short name= UEV-3 Alternative name(s): EV and lactate/malate dehydrogenase domain-containing protein
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.