



Recombinant Danio rerio UPF0406 protein C16orf57 homolog (zgc:91896)

Product Code	CSB-BP735289DIL
Abbreviation	zgc:91896
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DEF6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MIVNYSSSSS EEESGSSSSP SGKRQKLDTE TSEALDHGSA QRKVCKSSHL TPRLPLPESV KEMFRDSEEL WTDKSEEHGG RLRSFQHERG NWATYVFFPY DPEEAFLEVL NKMMAAAEAH DIPLTVSEEF HLSLSKTVVL RHHWIQPFVQ SIRTSLTHFQ KFYCVAYKLIK VYSNAEKTRT FLGMEVSTGT PHLLELSKIV DETMKEFNLS TFYEDPSFHI SLAWCVGDQT ERLKKACLLE LQGLIDAHED GPFHARLNCN ELRCKTGKNKV FVFPLQ
Source	Baculovirus
Target Names	usb1
Protein Names	Recommended name: UPF0406 protein C16orf57 homolog
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.