



Recombinant Rat Abhydrolase domain-containing protein 14B (Abhd14b)

Product Code	CSB-BP735304RA
Abbreviation	Abhd14b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DGG1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MANVELSEGT IQVRGQSLFF REARPGNGQA VRFSVLLHGH IRFSSETWQN LGTLHRLAEA GYRAVAIDL GLGRSKEAAA PAPIGELAPG SFLAAVVDAL ELGSLVVISP SLSGMYALPF LVAPESQLRG YVPVAPICTD KINAADYARV KTPTLIVYGD QDPMGSSSFQ HLKQLPNHRV LVMEGAGHPC YLDKPDEWHT GLLDLFLQELA
Source	Baculovirus
Target Names	Abhd14b
Protein Names	Recommended name: Abhydrolase domain-containing protein 14B EC= 3.-.-
Expression Region	1-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.