



Recombinant Danio rerio NF-kappa-B inhibitor-interacting Ras-like protein 2 (nkiras2)

Product Code	CSB-BP735305DIL
Abbreviation	nkiras2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DGL2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MGKSCKVVVC GQGSVGKTAV LEQLLYANHV VGSETMETLE DIYIGSVETD RGTREQVRFY DTRGLRDGQE FPRHYFTFAD GFVLVYSIDS RESFKRVEAL KKEIDRCRDK KEVTIVVVG N KLDLQDQRRV DSSAAQQWAR QEKVRLWEIS VTDRRTLIEP FVHLASKMTQ PQSKSTFPLS RNKNKTSGSL DS
Source	Baculovirus
Target Names	nkiras2
Protein Names	Recommended name: NF-kappa-B inhibitor-interacting Ras-like protein 2 Alternative name(s): I-kappa-B-interacting Ras-like protein 2 Short name= Kappa B-Ras protein 2 Short name= KappaB-Ras2
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.