



Recombinant Mouse Exocyst complex component 3-like protein 4 (Exoc3l4)

Product Code	CSB-BP735315MO
Abbreviation	Exoc3l4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DIA2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPLPQTGAPG PEVKSPREPR KSQTLPVTTW KSNSMKEQSV HHGGSLRPSL GMLKQTLFRT SLRTSTHKPK EDPGLFRRSS RFLFRSLRRA IDEGLTAGHP QGPAVPEKPS KVTDGVSQRQA ATGTEAEDLE PQAESKSVAD LITERQLVKA FEQLRYLETQ LVADKTSRTF TQDPTAYARR AMDLCLHYDG MAAEIGAIVR EALSSEGVD R DALAELAQVV HLEEEAHQTS QAEGDFLSTP RHWRMHWEDA VRLSAQESVQ QAGAKVIPGA AEGSSDLAQL LAELGGVVRH DLQKVRLEMQ PAYEATDFPV WETYLRAFHS AVAQRQLQELA RDARGCEQLY VLLDWAANVY GSPDFLGAPD LALPTEPLPP LLEPALWARL ESDYTSFLET KITSCFDSIL QLEQNRWEAD EDREVLQGLY HAPLSIDVHM LVAEHVKAAG AISAELEATT LQICARALCL FVPRFEKAFL ASKAVSEWYL GAYINACVEL RTSLLARFPG TIKELEKPLV AATNSFQKHL LQIVQQDMQP LFKVLYTKSW LTQDTLRPLM DKVVDFAHHL EHVTPPLAQE TLQEVRHFVV REYLGQVLRP HERFSGQDRL KGSNKMNDA QAISNTFQGL GSEAKWLDQA ILSVAEILGE TYKDDIRRH L ETLIRSYDPDI RRDHILAILA LRRLGRRRNQ NLLQHTQDLL RAAHETRLPS HHVLFEEIEV PTSVDVLITC I
Source	Baculovirus
Target Names	Exoc3l4
Protein Names	Recommended name: Exocyst complex component 3-like protein 4
Expression Region	1-721
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.