



Recombinant Mouse Exocyst complex component 3-like protein 4 (Exoc3l4)

| | |
|--------------------------|--|
| Product Code | CSB-YP735315MO |
| Abbreviation | Exoc3l4 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6DIA2 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MPLPQTGAPG PEVKSPREPR KSQTLPVTTW KSNSMKEQSV HHGGSLRPSL GMLKQTLFRT SLRTSTHKPK EDPGLFRRSS RFLFRSLRRA IDEGLTAGHP QGPAVPEKPS KVTDGVSRQA ATGTEAEDLE PQAESKSVAD LITERQLVKA FEQLRYLETQ LVADKTSRTF TQDPTAYARR AMDLCLHYDG MAAEIGAIVR EALSSEGVD R DALAELAQVV HLEEEAHQTS QAEGDFLSTP RHWRMHWEDA VRLSAQESVQ QAGAKVIPGA AEGSSDLAQL LAELGGVVRH DLQKVRLEMQ PAYEATDFPV WETYLRAFHS AVAQRQLQELA RDARGCEQLY VLLDWAANVY GSPDFLGAPD LALPTEPLPP LLEPALWARL ESDYTSFLET KITSCFDSIL QLEQNRWEAD EDREVLQGLY HAPLSIDVHM LVAEHVKAAG AISAELEATT LQICARALCL FVPRFEKAFL ASKAVSEWYL GAYINACVEL RTSLLARFPG TIKELEKPLV AATNSFQKHL LQIVQQDMQP LFKVLYTKSW LTQDTRLPLM DKVVDFAHHL EHVTPPLAQE TLQEVRHFVV REYLGQVLRP HERFSGQDRL KGSNKMNDA QAISNTFQGL GSEAKWLDQA ILSVAEILGE TYKDDIRRH L ETLIRSY PDI RRDHILAILA LRRLGRRRNQ NLLQHTQDLL RAAHETRLPS HHVLFEEIEV PTSVDVLITC I |
| Source | Yeast |
| Target Names | Exoc3l4 |
| Protein Names | Recommended name: Exocyst complex component 3-like protein 4 |
| Expression Region | 1-721 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final |



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.