



# Recombinant Macaca mulatta Renin (REN)

<b>Product Code</b>	CSB-MP735332MOW
<b>Abbreviation</b>	REN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6DLW5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Macaca mulatta (Rhesus macaque)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	LALG NTTSSVILTN YMDTQYYGEI GIGTPPQTFK VVFDTGSSNV WVPSSKCSRL YTACVYHKLF DASDSSSYKH NGTELTLRYS TGTVSGFLSQ DIITVGGITV TQMFGEVTEM PALPFMLAEF DGVVGMGFIE QAIGRVPIF DNILSQGVLK EDVFSFYNR DSENAQSLGG QIVLGGSDPQ HYEGNFHYIN LIKTGVWQIP MKGVS VGSST LLCEDGCLAL VDTGASYISG STSSIEKLME ALGAKKRLFD YVVKCNEGPT LPDISFHLGG KEYTLTSADY VFQESYSSKK LCTLAIHAMD IPPPTGPTWA LGATFIRKIFY TEFDRRNNRI GFALAH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	REN
<b>Protein Names</b>	Recommended name: Renin EC= 3.4.23.15 Alternative name(s): Angiotensinogenase
<b>Expression Region</b>	67-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Renin catalyzes the first step in the activation pathway of angiotensinogen--a cascade that can result in aldosterone release, vasoconstriction, and increase in blood pressure. Renin, an aspartyl protease, cleaves angiotensinogen to form angiotensin I, which is converted to angiotensin II by angiotensin I converting enzyme, an important regulator of blood pressure and electrolyte balance. Transcript variants that encode different protein isoforms and that arise from alternative splicing and the use of alternative promoters have been described, but their full-length nature has not been determined. Mutations in this gene have been shown to cause familial hyperproreninemia.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.