



Recombinant *Phytolacca heterotepala* Heterotepalin-4 (RIP1)

Product Code	CSB-EP735359PGAO-B
Abbreviation	RIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6EH50
Product Type	Recombinant Protein
Immunogen Species	<i>Phytolacca heterotepala</i> (Mexican pokeweed)
Purity	≥85% (SDS-PAGE)
Sequence	VNTIIYNV GSTTISKYAT FLDDLNRNEAK DPNLKCYGIP MLPNTNSNPK YVLVELQGSN KKTITLMLRR NNLYVMGYSD PFDTSKCRYH IFNDISGTER QDVETTLCPN SNSRVSKNIN YDSRYPTLES KVGVKRSRSQV QLGIQILDSD IGKISGVTSF SEKTEAEFLL VAIQISEAAR FKYIENQVKT NFNRAFNPNP KVLNLEETWG KISTAIHDAK NGVLPKPLEL VDASGAKWIV LRVDEIKPDV ALLNYVSGSC QTT
Source	E.coli
Target Names	RIP1
Protein Names	Recommended name: Heterotepalin-4 EC= 3.2.2.22 Alternative name(s): Ribosome-inactivating protein rRNA N-glycosidase
Expression Region	23-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.