



# Recombinant Anaplasma marginale tRNA-specific 2-thiouridylase mnmA (mnmA)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP735430AAAC                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Abbreviation</b>      | mnmA                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                              |
| <b>Uniprot No.</b>       | Q5PBB1                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Anaplasma marginale (strain St. Maries)                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Sequence</b>          | MQFIPSPAMV FDGSLNLDPL VSGKPPEETT VVAMSGGVD SSVVAALLHE<br>RGYKVIGATM QLHSSSPASG AKSCCGSVDI YDAKRVASTL GFPHYVLDYE<br>EVFRREVIDD FINSYKRGET PIPCVKCNQT VKFRDMLKAA RAIGGDVVAT<br>GHYVRRVEIA GEQQILRGKD PQKDQSYFLF SVTSEQLKFL RFPLGDLAKS<br>NVRLLAQQLN LEVADKPDSQ DICFVPENSY REVLRLNDPA SVKKKGKIVHV<br>DGRLLGEHDG ISNFTVGQRR GLNISAPYPL YVVRLDAAQN TVVVG PQSAL<br>MKRALFVKNL NWLPDYNIPK RGLAVDARLR SSGATVRATI TSDGDGYGTV<br>VLEEDCVVSP GQACVLYDKE RLLGGGWIYN KLCHEGA |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Target Names</b>      | mnmA                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | Recommended name: tRNA-specific 2-thiouridylase mnmA EC= 2.8.1.-                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Expression Region</b> | 1-387                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                             |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                     |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                              |