



Recombinant Salmonella paratyphi A GTP cyclohydrolase-2 (ribA)

Product Code	CSB-MP735443SBM
Abbreviation	ribA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PCY8
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MQLKRVAEAK LPTPLGDFLM VGFEELATGH DHAALVFGDI SGKTPVLARV HSECLTGDAL FSLRCDGCFQ LEAALTHIAE EGRGILIYHR QEGRNIGLLN KIRAYALQDQ GYDTVEANHQ LGFAADERDF TLCADMFKLL GVDEVRLLTN NPKKVEILTE AGINIVERVP LIVGRNPNE HYLDTKAAKM GHLLSK
Source	Mammalian cell
Target Names	ribA
Protein Names	Recommended name: GTP cyclohydrolase-2 EC= 3.5.4.25 Alternative name(s): GTP cyclohydrolase II
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.