



Recombinant Salmonella paratyphi A Ribosomal protein S12 methylthiotransferase RimO (rimO)

Product Code	CSB-MP735473SBM
Abbreviation	rimO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PGP7
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MSNVTHQPKI GFVSLGCPKN LVDSEIRLTE LRTEGYDVVP RYDDADMVIV NTCGFIDSAV QESLEAIGEA LNENGKVIVT GCLGAKEDQI REVHPKVLEI TGPHSYEQVL QHVHHYVPKP KHNPFLSLVP EQGVKLTPRH YAYLKISEGC NHRCTFCIIP SMRGDLVSRP IGDVLSEAKR LVDAGVKEIL VISQDTSAYG VDVKHRTGFH NGEPVKTSMV SLCEQLSKLG VWTRLHYVYP YPHVDDVIPL MAEGKILPYL DIPLQHASPR ILKLMKRP GS VDRQLARIKQ WREICPELTL RSTFIVGFPG ETEEDFQMLL DFLKEARLDR VGCFKYSPE GAGANELPDQ VPEEVKEERW NRFMLQLQQI SAERLQEKVG REILVIVDEV DEEGAIGRSM ADAPEIDGAV YLNGETNVKL GDIVRVKVEN ADEYDLWGSR V
Source	Mammalian cell
Target Names	rimO
Protein Names	Recommended name: Ribosomal protein S12 methylthiotransferase RimO Short name= S12 MTTase Short name= S12 methylthiotransferase EC= 2.-.-.- Alternative name(s): Ribosome maturation factor RimO
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.