



# Recombinant Salmonella paratyphi A Ribonuclease BN (rbn)

<b>Product Code</b>	CSB-EP735537SBM-B
<b>Abbreviation</b>	rbn
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5PN69
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MELIFLG TSA GVPTRSRNVT AILLHLQHPT QPGVWLFDCG EGTQHQMLNT AFHPGKLERI FISHLHGDHL FGLPGLLCSR SMAGNPHPLT VYGPQGVREF IATTLRLSGS WTD FPLQIEE ISAGDILDDG LRVKTA FRLE HPLECYGYRV VEHDKPGALN ARALKAAGVT PGPLFQALKA GKT VTLADGR QINGADYLAP AVAGKSVAIF GDTAPCEAAL ALAQQVDVMV HETTL DASME EKANARGHSS TRQTATLARE AAVGRLIMTH ISSRYDDKGC QRLLAECRAI FPATELAYDF SVFPV
<b>Source</b>	E.coli
<b>Target Names</b>	rbn
<b>Protein Names</b>	Recommended name: Ribonuclease BN Short name= RNase BN EC= 3.1.-.- Alternative name(s): Ribonuclease Z homolog Short name= RNase Z homolog
<b>Expression Region</b>	1-305
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.