



Recombinant *Oryza sativa* subsp. *japonica* Tubby-like F-box protein 1 (TULP1)

Product Code	CSB-BP735577OFG
Abbreviation	TULP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5QM27
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MPRETPPPPP PEGGEVHEVV EGEDGQAEDQ EERWARLLPE LMSEVRRVE ASGGERWPAR KDVVSCACVC RRWRDAAVAV VRPPAESGKI TFPSSLKQPG PREFPMQCFI KRKKNSTFY LYLGLTNATV DKGKFLMAAR RFRRGPHTY IVSLDADDLS QGSNAYMGKL RSDFWGTNFK IYDSKPPYDG AKASSSRSSR RFGSRRISPQ VSAGNYEVGQ VSYKYNLLKS RGPRRMYCAL ECPSTQETWE NCLKTKFRKP TGNTVLRNKA PRWHEHLQCW CLNFHGRVTV ASVKNFQLVA AADPNDPASS KDEETVLLQF GKVDDNIFTM DYRQPLSAFQ AFAISLSSFG TKLACE
Source	Baculovirus
Target Names	TULP1
Protein Names	Recommended name: Tubby-like F-box protein 1 Short name= OsTLP1 Alternative name(s): Tubby-like F-box protein 9 Short name= OsTLP9
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.