



Recombinant Pongo abelii Septin-7 (41524)

Product Code	CSB-BP735643PYX
Abbreviation	41524
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R481
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	SVSARSAAA EERSVNSSTM VAQQKNLEGY VGFANLPNQV YRKSVKRGFE FTLMVVGESG LGKSTLINS L FLTDLYSPEY PGPSHRIKKT VQVEQSKVLI KEGGVQLLLT IVDTPGFGDA VDNSNCWQPV IDYIDSKFED YLNAESRVNR RQMPDNRVQC CLYFIAPSGH GLKPLDIEFM KRLHEKVNII PLIAKADTLT PEECQQFKKQ IMKEIQEHKI KIYEFPETDD EEENKLVKKI KDRLPLAVVG SNTIIEVNGK RVRGRQYPWG VAEVENGEHC DFTILRNMLI RTHMQDLKDV TNNVHYENYR SRKLAAVTYN GVDNNKNKGQ LTKSPLAQME EERREHVAKM KKMEMEMEQV FEMKVKEKVQ KDKDSEAE LQ RRHEQMKKNL EAQHKELEEK RRQFEDEKAN WEAQQRILEQ QNSSRTLEKN KKKGKIF
Source	Baculovirus
Target Names	SEPT7
Protein Names	Recommended name: Septin-7
Expression Region	2-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.