



Recombinant Pongo abelii Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma (PIP4K2C)

Product Code	CSB-MP735644PYX
Abbreviation	PIP4K2C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R488
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	ASSSVPPAT VSAATTGPGP GFGFASKTKK KHVQEQEVKV FRAADPLVGV FLWGVVAHSIN ELSQVPPPVM LLPDDFKASS KIKVNNHLFH RENLPSHFKE KEYCPQVFRN LRDRFGIDDQ DYLVSLTRNP PSESEGSDGR FLICYDRTL IKEVSSEDIA DMHSNLSNYH QYIVKCHGNT LLPQFLGMYR VSVDNEDSYM LVMRNMFSHR LPVHRKYDLK GSVSREASD KEVKELPTL KGMDFLNKNQ KVYIGEDEKK IFLEKLKRDV EFLVQLKIMD YSLLLGIHDI IRGSEPEEEG PVREDESEVD GDCSLTGPPA LVGSYGTSPE GIGGYIHSR PLGPGEFESF IDVHAIRSAE GAPQKEVYFM GLIDILTQYD AKKKAHAHA AK TVKHGAGAEI STVHPEQYAK RFLDFITNIF A
Source	Mammalian cell
Target Names	PIP4K2C
Protein Names	Recommended name: Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma EC= 2.7.1.149 Alternative name(s): Phosphatidylinositol 5-phosphate 4-kinase type II gamma Short name= PI(5)P 4-kinase type II gamma Short name= PIP4KII-
Expression Region	2-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.