



Recombinant Pongo abelii MOB kinase activator 3A (MOB3A)

Product Code	CSB-EP735683PYX-B
Abbreviation	MOB3A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R5Z0
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MSNPFLKQVF NKDKTFRPKR KFEPGTQRFE LHKKAQASLN AGLDLRLAVQ LPPGEDLNDW VTVHVVDFFN RVNLIYGTIS DGCTEQSCPV MSGGPKYEYR WQDEHKFRKP TALSAPRYMD LLMDWIEAQI NNEDLFPTNV GTPFPKNFLQ AVRKILSRLE RVFVHVYIHH FDRIAQMGSE AHVNTCYKHF YYFVKEFGLI DTKELEPLKE MTARMCH
Source	E.coli
Target Names	MOB3A
Protein Names	Recommended name: MOB kinase activator 3A Alternative name(s): Mob1 homolog 2A Mps one binder kinase activator-like 2A
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.