



Recombinant Pongo abelii COMM domain-containing protein 2 (COMMD2)

Product Code	CSB-MP735685PYX
Abbreviation	COMMD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R610
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MLLELSEEHK EHLAFLPQVD SSVVAEFGRI AVEFLRRGAN PKIYEGAARK LNVSSDTVQH GVEGLTYLLT ESSKLMISEL DFQDSVFLG FSEELNLLLL QLYLDNRKEI RMILSELAPS LPSYHNLEWR LDVQLASRSL RQGIKPAVTI KLHLNQNQGDH NTKILQTDPA TLLHLVQQLE QALEEMKTNH CRRVVRNIK
Source	Mammalian cell
Target Names	COMMD2
Protein Names	Recommended name: COMM domain-containing protein 2
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.