



Recombinant Pongo abelii Platelet-activating factor acetylhydrolase IB subunit gamma (PAFAH1B3)

Product Code	CSB-BP735701PYX
Abbreviation	PAFAH1B3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R6X1
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	SGEENPASK PTPVQDVQGD GRWMSLHHRF VADSKDKEPE VVFIGDSLVSQ LMHQCEIWRE LFSPLHALNF GIGGDGTQHV LWRENGELE HIRPKIVVWV VGTNNHGHTA EQVTGGIKAI VQLVNERQPQ ARVVVLDLLP RGQHPNPLRE KNQRVNELVR AALAGHPRAH FLDADPGFVH SDGTISHHDM YDYLHLSRLG YAPVCRALHS LLLRLLAQDQ GQGAPLLDPA P
Source	Baculovirus
Target Names	PAFAH1B3
Protein Names	Recommended name: Platelet-activating factor acetylhydrolase IB subunit gamma EC= 3.1.1.47 Alternative name(s): PAF acetylhydrolase 29 kDa subunit Short name= PAF-AH 29 kDa subunit PAF-AH subunit gamma Short name= PAFAH su
Expression Region	2-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.