



Recombinant Pongo abelii MOB kinase activator 1A (MOB1A)

Product Code	CSB-BP735761PYX
Abbreviation	MOB1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RAE0
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	SFLFSSRSS KTFKPKKNIP EGSHQYELLK HAEATLGSGN LRQAVMLPEG EDLNEWIAVN TVDFFNQINM LYGTITEFCT EASCPVMSAG PRYEYHWADG TNIKKPIKCS APKYIDYLMT WVQDQLDDET LFPSKIGVPF PKNFMSVAKT ILKRLFRVYA HIYHQHFDSV MQLQEEAHLN TSFKHFIFFV QEFNLIDRRE LAPLQELIEK LGSKDR
Source	Baculovirus
Target Names	MOB1A
Protein Names	Recommended name: MOB kinase activator 1A Alternative name(s): Mob1 homolog 1B Mps one binder kinase activator-like 1B
Expression Region	2-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.