



Recombinant Pongo abelii Axin interactor, dorsalization-associated protein (AIDA)

Product Code	CSB-EP735770PYX
Abbreviation	AIDA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RAV3
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MSEVTRSLQ RWGASFRRGA DFDSWGQLVE AIDEYQILAR HLQKEAQAQH NNSEFTEEK KTIGKIATCL ELRSAALQST QSQEETFKLED LKKLEPILKN ILTYNKEFPF DVQPVPLRRI LAPGEEENLE FEEDEEEGGA GAGSPDSFPA RVPGTLLPRL PSEPGMTLLT IRIEKIGLKD AGQCIDPYIT VSVKDLNGID LTPVQDTPVA SRKEDTYVHF NVDIELQKHV EKLTKGAAIF FEFKHYKPKK RFTSTKCFAP MEMDEIKPGP IVIELYKKPT DFKRKKLQLL TKKPLYLHLH QTLHKE
Source	E.coli
Target Names	AIDA
Protein Names	Recommended name: Axin interactor, dorsalization-associated protein Alternative name(s): Axin interaction partner and dorsalization antagonist
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.