



Recombinant Pongo abelii Coiled-coil domain-containing protein 86 (CCDC86)

Product Code	CSB-BP735776PYX
Abbreviation	CCDC86
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RB69
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	MDTPLRRSRR LGGLRPECPE SLTSVSRTRR ALVELESNPE ETREPGSPPS VQRAGLGSPE RPPKTSPGSP RLQQGSGLES PQQQPEPGAG SPQRQQDLHL ESPQRQPEYS PESPRCQPKP SEEVPKCSQD QGVLASELAQ SKEELTPGAP QHQLPPVPGS PEPYPGQQAP GPEPSQPLLE LTPRAPGSPR GQHEPSKPPP AGETVTGGFG AKKRKGSSSL APASKKLNKE ELPVIPKGKP KSGRVWKDRS KKRFSQMLQD KPLRTSWQRK MKERQERKLA KDFARHLEEE KERRRQEKKQ RRAENLKRRL ENERKAEVVQ VIRNPAKLKR AKKKQLRSIE KRDTLALLQK QPPQRPAAKI
Source	Baculovirus
Target Names	CCDC86
Protein Names	Recommended name: Coiled-coil domain-containing protein 86
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.