



Recombinant Pongo abelii EH domain-containing protein 1 (EHD1), partial

Product Code	CSB-MP735788PYX
Abbreviation	EHD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RBP4
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	MFSWVSKDAR RKKEPELFQT VAEGLRQLYA QKLLPLEEHY RFHEFHSPAL EDADFDNKPM VLLVGQYSTG KTFIRHLIE QDFPGMRIGP EPTTDSFIIV MHGPTGVVP GNALVVDPRR PFRKLNAFGN AFLNRFMCAQ LPNPVLDSIS IIDTPGILSG EKQRISRGYD FAAVLEWFAE RVDRIILLFD AHKLDISDEF SEVIKALKNH EDKIRVVLNK ADQIETQQLM RYVYALMWSL GKIINTPEVV RVYIGSFWSH PLLIPDNRKL FEAEEQDLFK DIQSLPRNAA LRKLNLIKR ARLAKVHAYI ISSLKKEMPN VFGKESKKKE LVNNLGEIYQ KIEREHQISP GDFPSLRKMQ ELLQTQDFSK FQALKPKLLD TVDDMLANDI ARLMVMVRQE ESLMPSQVVK GGAFDGTMNG PFGHGYGEGA GEGIDDVEWV VGKDKPTYDE IFYTLSPVNG KITGANAKKE MVKSKLPNTV LGKIWKLADV DKDGLLDDEE FALANHLIKV KLEGHELPAD LPPHLVPPSK RRHE
Source	Mammalian cell
Target Names	EHD1
Protein Names	Recommended name: EH domain-containing protein 1
Expression Region	1-534
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.