



# Recombinant Pongo abelii EH domain-containing protein 1 (EHD1), partial

<b>Product Code</b>	CSB-YP735788PYX
<b>Abbreviation</b>	EHD1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5RBP4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFSWVSKDAR RKKEPELFQT VAEGLRQLYA QKLLPLEEHY RFHEFHSPAL EDADFDNKPM VLLVGQYSTG KTFIRHLIE QDFPGMRIGP EPTTDSFIIV MHGPTGVVP GNALVVDPRR PFRKLNAFGN AFLNRFMCAQ LPNPVLDSIS IIDTPGILSG EKQRISRGYD FAAVLEWFAE RVDRIILLFD AHKLDISDEF SEVIKALKNH EDKIRVVLNK ADQIETQQLM RYVYALMWSL GKIINTPEVV RVYIGSFWSH PLLIPDNRKL FEAEQDLFK DIQSLPRNAA LRKLNLIKR ARLAKVHAYI ISSLKKEMPN VFGKESKKKE LVNNLGEIYQ KIEREHQISP GDFPSLRKMQ ELLQTQDFSK FQALKPKLLD TVDDMLANDI ARLMVMVRQE ESLMPSQVVK GGAFDGTMNG PFGHGYGEGA GEGIDDVEWV VGKDKPTYDE IFYTLSPVNG KITGANAKKE MVKSKLPNTV LGKIWKLADV DKDGLLDDEE FALANHLIKV KLEGHELPAD LPPHLVPPSK RRHE
<b>Source</b>	Yeast
<b>Target Names</b>	EHD1
<b>Protein Names</b>	Recommended name: EH domain-containing protein 1
<b>Expression Region</b>	1-534
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.