



Recombinant Pongo abelii Neuronal calcium sensor 1 (NCS1)

Product Code	CSB-BP735799PYX
Abbreviation	NCS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RC90
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	GKSN SKLKP EVVEELTRKT YFTEKEVQQW YKGF IKDCPS GQLDAAGFQK IYKQFFPFGD PTKFATFVFN VFDENKDGRI GFSEFIQALS VTSRGTLDEK LRWAFKLYDL DNDGYITRNE MLDIVDAIQ MVGNTVELPE EENTPEKRVD RIFAMMDKNA DGKLT LQEFQ EGSKADPSIV QALSLYDGLV
Source	Baculovirus
Target Names	NCS1
Protein Names	Recommended name: Neuronal calcium sensor 1 Short name= NCS-1 Alternative name(s): Frequenin homolog
Expression Region	2-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.