



Recombinant Pongo abelii Probable RNA-binding protein EIF1AD (EIF1AD)

Product Code	CSB-MP735817PYX
Abbreviation	EIF1AD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RD29
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MSQATKRKHV VKEVLGEHIV PSDQQQIVRV LRTPGNNLHE VETAQGQRFL VSMPSKYRKN IWIKRGDFLI VDPIEEGEKV KAEISFVLCK DHVRS LQKEG FWPEAFSEVA EKHNNNRNRQ TQPELPAEPQ LSGEESSSED DSDLFVNTNR RQYRESEEEES EEEEEAA
Source	Mammalian cell
Target Names	EIF1AD
Protein Names	Recommended name: Probable RNA-binding protein EIF1AD Alternative name(s): Eukaryotic translation initiation factor 1A domain-containing protein
Expression Region	1-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.