



Recombinant Pongo abelii Probable methyltransferase BCDIN3D (BCDIN3D)

Product Code	CSB-BP735871PYX
Abbreviation	BCDIN3D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RFI3
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MAVPTELHGG SVKETAAEKE SRVLEPGAAP FGNFPHYSRF HPPEQRLRL PPELLRQLFP ESPENGPILG LDVGCNSGDL SVALYKHFLS LPDGKTCSDA SREFRLLCCD IDPVLVKRAE KECFPFDALT FITLDFMNQR TRKVLLSSFL SQFGRSVFDI GFCMSITMWI HLNHGDHGLW EFLARLSSLC RYLLVEPQPW KCYRAAARRL RKLGLHDFDH FHSLTIRGDM PNQIVQILTQ DHGMELICCF GNTSWDRSLL LFRAKQTIET HPIPESLIEK GKEKNRLSFQ KQ
Source	Baculovirus
Target Names	BCDIN3D
Protein Names	Recommended name: Probable methyltransferase BCDIN3D EC= 2.1.1.- Alternative name(s): BCDIN3 domain-containing protein
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.