



Recombinant Danio rerio Transcription factor AP-2-delta (tfap2d)

Product Code	CSB-EP735885DIL-B
Abbreviation	tfap2d
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RJ20
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSATFPGLVH DAEIRHDGSN SYRLMQLGCL ESVANSSVAY SSSSPLTYPA TGTEFASPYF PTNHQYTPLH HQSFHYEFQH SHPAVNPDAY SLNSLHHSQQ YYQQLHHGEP ADFINLHNAR ALKSSCLDEQ RRELGCLDAY RRHDLSIMSH GSQYGMHPEQ RLLPGAGLGL PPPGADDLQG SVDAQCGLVL NGQGGVIRRG GTCVVNPTDL FCSVPGRLSL LSSTSKYKVT IAEVKRRLSP PECLNASLLG GILRRAKSKN GGRCLREKLD RLGLNLPAGR RKAANVTLT SLVEGEALHL ARDFGYTCET EFPTKAVGEH LARQHTEPKE QNARKKMVLA TKQICKEFQD LLSQDRSPLG SSRPTPILD DLQIRHLTHFS LITHGFGTPA ICAALSTFQT ILSEMLNYLE KHSANKSTGT PDTNQINSNS DKTLRKTEAP TKDGKIEKTE
Source	E.coli
Target Names	tfap2d
Protein Names	Recommended name: Transcription factor AP-2-delta Short name= AP-2-delta
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.