



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 67 (At5g63260)

Product Code	CSB-EP735886DOA-B
Abbreviation	At5g63260
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RJC5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSKPEETSDP NPTGPDPSRS SSDEVTVTVA DRAPSDLNHV SEELSDQLRN VGLDDSAKEL SVPIVPPQGN VETDSRALFG SDQKEEEEGS EKRMVVYPV RPDSEDCSFY MRTGSCYKYG SCKFNHPVRR KLQIGRERVR ERDEDVENPK LMECKYYFRT GGCKYGESCR FSHMKEHNSP ASVPELNFLG LPIRPGEKEC PFYMRNGSCK FGSDCKFNHP DPTAIGGVDS PLYRGNNGGS FSPKAPSQAS STWSSTRHM NGTGTAPFIP SMFPHSRGVT PQASDWNGYQ ASSAYPPERS PLAPSSYQVN NSLAETSSFS QYQHQMVSVEE FPERPDQPEC TYYLKTGDCK FKYKCKYHHP KNRLPKQAAF SFNDKGLPLR PDQSMCTHYS RYGICKFGPA CRFDHSIPPT FSPSSSQTVE ARQVGANGNE DDSWH
Source	E.coli
Target Names	At5g63260
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 67 Short name= AtC3H67 Alternative name(s): Zinc finger CCCH domain-containing protein ZFN-like 5
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.