



# Recombinant Chicken Polycomb complex protein BMI-1 (BMI1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP735906CH  |
| <b>Abbreviation</b>      | BMI1  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q5SDR3  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Gallus gallus (Chicken)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MHRTTRIKIT ELNPHLMCVL CGGYFIDATT IIECLHSFCK TCIVRYLETS<br>KYCPICDVQV HKTRPLLNIR SDKTLQDIVY KLVPGLFKNE MKRRRDFYAA<br>HPSADAANGS NEDRGEVADE DKRIITDDEI ISLSIEFFDQ NRLERKGNKE<br>KEKSKEEVND KRYLRCPAAM TVMHLRKFLR SKMDIPNTFQ IDVMYEEEEPL<br>KDYYTLMDIA YIYTWRNGP LPLKYRVRPT CKRMKIGHQR EGLSNSGELE<br>SDSGSDKASS PAGGLPSTSS CLPSPSTPVQ SPHPQFPHIS STMNGTSSSP<br>GSNHQSSFTN RARKSSINGS SATSSG                    |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | BMI1  |
| <b>Protein Names</b>     | Recommended name: Polycomb complex protein BMI-1 Alternative name(s): Polycomb group RING finger protein 4  |
| <b>Expression Region</b> | 1-326   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |